

There is a global shortage of Tru Steel infusion sets which are no longer available.

Patients must switch to a new infusion set on Tandem Pump.

Recommendations for ease of use are in order below (easiest starting with vari soft):

Vari Soft: <https://diabetesexpress.ca/collections/tandem-t-lock-insulin-pump-supplies/products/t-lock-varisoft-tandem-infusion-set>

- a. Video: https://youtu.be/gPml_HjqFPk?si=wIOzQqcLT6pwmvcr

Auto Soft 30 (AS30): <https://diabetesexpress.ca/collections/tandem-t-lock-insulin-pump-supplies/products/t-lock-autosoft-30-tandem-infusion-set>

- a. https://youtu.be/sQJcab_5FCw?si=bRsE4FoJEzQfRSV

Autosoft 90 (AS90): <https://diabetesexpress.ca/collections/tandem-t-lock-insulin-pump-supplies/products/t-lock-autosoft-90-tandem-infusion-set-1>

- a. <https://youtu.be/y-EBQ3QgMaw?si=6KtWd64ZawQQ5KPy>